

Technical data:

PARAMETER:	AR-201	AR-201A	AR-202	AR-203	AR-204	AR-205	AR-208	AR-209
Frame:	aluminum	aluminum	aluminum	aluminum	steel	aluminum	aluminum	aluminum
Frame color:	silver	silver/gray	silver	silver	white	silver	silver	silver
Type/model	constant	constant	constant	constant	folding	folding	constant	constant
Seat width:	33 cm round	33 cm round	50 cm	50 cm	50 cm	50 cm	48 cm	41 cm
Seat depth:	33 cm round	33 cm round	30 cm	30 cm	30 cm	30 cm	33 cm	48 cm
Height adjustment	38-56 cm	37-52 cm	36-47cm	36-47 cm	52 cm	47-50.5 cm	36-53 cm	35-45 cm
Product weight:	1.5kg	1.9 kg	1.8 kg	3.0kg	3.2kg	2.4kg	2.7kg	2.1 kg
Maximum load:	110 kg	110 kg	110 kg	110 kg	110 kg	110 kg	110 kg	110 kg



**USER MANUAL
WARRANTY CARD**

SHOWER CHAIRS AND STOOLS

SYMBOLS:

Wyrób medyczny	Zapoznaj się z instrukcją	Producent	Wyrób spełnia wymagania UE	Ostrzeżenie	Nr dato / serii	Data produkcji m-c/rok	Nadaje się do recyklingu	Selektywna zbiórka odpadów

WARRANTY CARD
The card must be accompanied by proof of purchase

Sale date	Seller's stamp and signature
Purchase invoice number	
Product no. / LOT no	

1. ARMEDICAL Sp. z o. o. 42-520 Dąbrowa Górnicza, al. Zwycięstwa 75 hereby grants a 12-month warranty for the products covered by this document. The warranty period is counted from the date of purchase of the product by the User. The warranty is valid only in the territory of the Republic of Poland.
2. During the warranty period, ARMEDICAL Sp. z o. o. undertakes to make the necessary repairs in order to restore the product to proper functioning within 14 days from the date of receipt of the complained goods.
3. ARMEDICAL Sp. z o. o. reserves the right to replace the product with a defect-free one if it considers that the repair costs make it unprofitable.
4. If it is found that the product is not functioning properly, please immediately contact the point of purchase or send it to the manufacturer's address.
5. The warranty does not cover damage caused by mechanical factors. The warranty does not cover the subject of sale in terms of wear and tear elements during normal use, i.e. scratches, scratches on polyethylene and rubber elements (i.e. abrasion of rubber caps (5)).
6. The warranty does not apply if the chair/shower stool is used contrary to its intended use or in accordance with the operating instructions, or if the user exceeds the maximum permissible load specified by the manufacturer.
7. The warranty does not cover defects resulting from unauthorized repairs, modifications or design changes to the product made by the User or other unauthorized persons.
8. **The warranty is valid only with the proof of purchase and this warranty card.**
9. Complaints should be submitted only in the form of a written complaint protocol specifying the type of defects found and the circumstances in which the defect was discovered. In case of failure to comply with the conditions for filing complaints specified in this manual, the guarantor's liability is excluded. **The advertised product must be in proper hygienic condition.**



- AR-202 Shower chair without backrest
- AR-203 chair with backrest
- AR-208 U-cut shower chair
- AR-201A Round shower stool - swivel

- AR-201 Shower stool - round
- AR-204 - foldable Shower stool
- AR-205 - foldable Shower stool with U cutout
- AR-209 - stool with U cutout



ARMEDICAL Sp. z o.o.
al. Zwycięstwa 75, 42-520 Dąbrowa Górnicza, Polska
tel.fax: +48 32 261 71 36, +48 783 001 111
www.armedical.pl, e-mail: biuro@armedical.pl





USER MANUAL - SHOWER CHAIRS AND STOOLS

Ladies and Gentlemen,

Thank you for purchasing our product, we value the trust you have placed in us. We have made every effort to ensure that our equipment meets your expectations, we hope that it will significantly facilitate your daily toileting activities. Before use, please read this user manual and strictly follow the safety regulations. Please keep this manual, which also serves as a warranty card.

Intended use of the product:

Shower chairs and stools as auxiliary equipment are intended for disabled people with disorders of the musculoskeletal system, elderly people as well as people for whom standard use of a bathtub or shower cabin is a problem for health reasons. They facilitate and help with everyday toilet activities, e.g. at the sink, in the bathtub or in the shower.



Contraindications:

Contraindications to a sitting position, lack of trunk stabilization. Movement coordination and balance disorders, labyrinth problems, spinal instability. Intellectually limited people should not use the product without the supervision of a guardian.

Usage data:

Maximum patient weight ensuring safe use of the product -**110 kg**



Notes and warnings:

To avoid an accident, please:

1. Checking whether all tubes are straight and not broken (e.g. during transport), and whether all product components are included in the packaging.
2. Checking whether all polyethylene elements have no damage, e.g. cracks
3. Checking whether all chair/stool legs are set at the same height. Applicable models: (AR-201, AR-201A, AR-202, AR-203, AR-205, AR-208, AR-209)
4. Check whether all legs have anti-slip rubber pads and whether they are all correctly installed.



To avoid injuries:

1. All polyethylene elements must be solidly and permanently attached (using the included screws) to the legs and backrest (applies to AR-203, AR-208 models)
2. It is necessary to check the technical condition and stability of the product each time before use (checking all fastening elements, product frame tubes and rubber caps)
3. All legs of the shower chair/stool must be set at the same height. Applicable models: (AR-201, AR-201A, AR-202, AR-203, AR-205, AR-208, AR-209)
4. Use only on a flat, hard surface that prevents the chair from moving when sitting down
5. Be especially careful on wet or slippery surfaces despite the non-slip caps
6. It is forbidden to sit down abruptly on a chair/stool (e.g. jump on it)
7. It is prohibited to move the chair/stool with the user
8. It is forbidden to stand with shoes on the seat or any element of the product.
9. Do not overload the product or exceed the maximum load specified by the Manufacturer.
10. Before using the product, it is recommended that slippery surfaces be covered with anti-slip mats (e.g. shower cabin).

Maintenance method:

Each time after toileting, the product should be thoroughly rinsed with clean water (do not leave any soap, shampoo, bath gel, etc. on the product) and then dried thoroughly with a delicate cloth. It is prohibited to use strong detergents, gasoline or solvents for cleaning. Storage in rooms with high humidity, e.g. (basement), combined with lack of thorough maintenance of the product (hygiene, drying) will result in loss of functional and aesthetic parameters of the chair/stool and thus loss of warranty. Do not use abrasive products! Anti-slip rubber elements (leg caps) (5) are wearing parts and must be regularly replaced with new ones.

Troubleshooting:

If you encounter problems, please review this chapter. If this does not help, please contact the sales or service point. Repairing the product yourself is FORBIDDEN!!! This may result in an accident.

I. The shower chair/stool is not stable.

- a. Are all chair/stool legs at the same height? If not, set them correctly.
- b. Do all legs have rubber caps (5)? If not, put them on.
- c. Are all rubber caps (5) installed correctly? If not, press them in firmly.
- d. Does the product have polyethylene elements firmly screwed to the frame (1) (with included screws)? If not, retighten all screws firmly. All product elements should be screwed/folded in such a way that, when completed, they constitute an exact reproduction of the product photo on page 1 (applies to all models).
- e. Despite firmly tightening all mounting screws, the product is still not stable? If not, please contact your point of purchase. DO NOT USE A SHOWER CHAIR/STOOL!!!

Incidents:

Everyoneserious an incident related to the device should be reported to the manufacturer, the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

Utilization:

A product that is no longer suitable for further use should be disposed of in accordance with local environmental regulations regarding product disposal.

Cessation of use of the product:

In case of visible mechanical damage or material wear, stop using the product.

COMPONENTS/UTILITY ELEMENTS



1. Aluminum frame
(Steel frame: applies to model AR-204)
2. Polyethylene seat
(shape depends on the model)
3. Polyethylene backrest
(applies to AR-203, AR-208 models)
4. Leg height adjustment
(applies to models AR-201, AR-201A, AR-202, AR-203, AR-205, AR-208, AR-209)
5. Anti-slip rubber caps

Preparation for use

1. Remove the shower chair/stool from the cardboard or foil packaging. The product removed from the packaging requires assembly before first use.
2. First, check the technical condition of the product in accordance with pointNotes and warnings.
3. Place the seat (2) on the floor or hard surface so that the bottom part is visible. Grab the elements of the aluminum frame (1) (steel frame for the AR-204 model) and place them with the rubber caps upwards (5) in such a way that the holes in the frame and seat elements coincide with each other. Screw all included screws into the existing holes, using your fingers - if using plastic butterfly screws. For the AR-203, AR-208 model: after assembling the seat and aluminum frame (legs), screw the backrest (3). Proceed as when assembling the seat and aluminum frame.
ATTENTION! All screws must be tightened tightly to ensure safe use of the product.
4. Adjust the height of the legs of the shower chair/stool (applies to models: AR-201, AR-201A, AR-202, AR-203, AR-205, AR-208, AR-209) to the user's height, for this purpose press the latch head (4) located in the leg to release the lock. Set the required height by extending or sliding the tube. The latch head must engage in the hole at the desired height! Repeat the operation for all legs. The latch head on all legs must be at the same height.

Alternatively, please consult your doctor, rehabilitator or point of sale for detailed instructions on appropriate use.

End of use:

1. Assemble the shower stool for models AR-204, AR-205

Cleaning and storage:

1. Place the shower chair/stool in a place where it will not be exposed to mechanical damage.
2. In case of dirt, follow the "Maintenance method" section of this manual.

THANK YOU FOR CHOOSING AND PURCHASING OUR PRODUCT